

More for the lab

Rolling cabinets for Hettich centrifuges free up space on the laboratory bench and give you greater flexibility.



Rolling Cabinets


All rolling cabinets are of sturdy design and feature castors and two locks. Any vibrations of the centrifuge are safely compensated for.

Rolling cabinets with Cat. No. 4612-A and 4613-A are suited for the models:

- ROTOFIX 32 A
- UNIVERSAL 320/320 R
- ROTINA 380/380 R
- ROTINA 420
- ROTANTA 460
- ROTOFIX 46/46 H

Rolling cabinets with Cat. No. 4614-K and 4615-A are suited for the models:

- ROTINA 420 R
- ROTANTA 460 R



Low and narrow,
with one drawer.
Depth 650 mm

Cat. No. 4612-A



Low and wide,
with one drawer.
Depth 650 mm

Cat. No. 4614-K



High and narrow,
with two drawers.
Depth 650 mm

Cat. No. 4613-A



High and wide,
with two drawers.
Depth 650 mm

Cat. No. 4615-A