

CE

# Coronavirus COVID-19



# Lambra

SPECIALISTS IN CLINICAL LABORATORY PRODUCTS SINCE 1989

EN

[www.Lambra.es](http://www.Lambra.es) | +34 91 798 56 43 | [Lambra@Lambra.com](mailto:Lambra@Lambra.com)


# Who are we?

*Lambra is a national company specialized in the commercialization of medical products since 1989, with over 30 years of experience in modernizing and improving the quality of our products and services. After many years of distributing hospital products for diagnosis, treatment and analysis, we continue to develop and create new concepts and ideas for our customers, guaranteeing the most complete systems for the identification and monitoring of your patients.*



# New product: Quick test COVID-19


*The quick test COVID-19 IgG/IgM (blood/serum/plasma) is an immunochromatographic test in a solid phase for the fast, qualitative and differential detection of antibodies IgG and IgM against the new coronavirus 2019 in whole blood, serum or human plasma.*






# Benefits and Functions:

- Results in 2-10 minutes.
- Facilitates patient treatment decisions quickly.
- Simple procedure + time saving.
- Small samples, 5 $\mu$ L of serum/plasma and 10 $\mu$ L of whole blood.
- Necessary reagents are provided without the need for medical equipment.
- High sensitivity and specificity.



# Process of detection and interpretation



# Lambra



+34 91 798 56 43



Travesía Igarsa,3  
28860  
Paracuellos de Jarama,  
Madrid



[lambra@lambra.com](mailto:lambra@lambra.com)



[Pedidos@Lambra.es](mailto:Pedidos@Lambra.es)



[www.Lambra.es](http://www.Lambra.es)

# Contact Us

