

Necklite™

Cervical Stabilization Device - Adult

- Cervical stabilization for extrication and rehabilitation
- Large tracheal opening allows for optimal visualization
- Helps reduce lateral movement
- Foldable and moldable before and after application
- One size fits all adults
- Single patient use
- X-ray translucent
- Fits in any kit
- Latex free
- Non-rigid



One Size Fits All Adults • No Settings • Device Molds to Patient

Introducing the Necklite™ - a moldable cervical stabilization device, comprised of soft aluminum housed within a high density PU foam. These characteristics allow the device to be folded or rolled for compact storage in backpacks, kits, and narrow pockets. A series of small, moldable tabs positioned under the chin conform to each patient's unique features. Once applied and adjusted, the device creates a rigid mold, securely stabilizing the neck.

Product Highlights

Light, compact design allows for virtually limitless storage options and convenient transport.

Fast, effective, and easy to use, even in demanding environments.

X-ray translucent design allows for cervical evaluation without removal of the collar.



The Necklite provides cervical stabilization in a low profile, compact package.

- Can safely be applied to patients found in awkward positions
- Molds to accommodate patient size
- Velcro® strap securing system
- Round front opening accommodates tracheotomy procedure
- Elliptical view port on back of device
- Rigid, X-ray translucent aluminum core
- Comfortable high density foam coating

Specifications

Dimensions 21.5" (550 mm) - Length
7" (180 mm) - Width
.2" (5 mm) - Thickness

Weight 3.7 oz (105 g)


Packaging Flat packed for easy shipping and storage

Sizes For all adult sizes

X-ray Translucent Yes

Available Colors Tactical / Military:
Dark Gray (both sides)

EMS / Responder:
Gray inside,
Yellow outside



PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com
www.ps-med.com

