



PRESSURE INFUSION CUFFS

Dispo Infusor 500 ml / 1000 ml / 3000 ml

Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes

Features

Transparent Netting

- Allows clear view on the fluid bag and the fluid level

Pocket Design with hook fastener/hook*

- Fluid bag can be replaced quickly and easily

Pressure Gauge

- Colour-coded with overpressure valve
- 360° visible

Hand Inflator

- Large volume allows quick inflation

Three-Way stopcock

- For pressure control
- For quick deflation

For single use

- Reduces the potential for cross-contamination

*With hook fastener: Dispo Infusor 500 ml and 1000 ml
With hook: Dispo Infusor 3000 ml



Dispo Infusor with female Luer Lock

Dispo Infusor with hand inflator and three-way stopcock

Dispo Infusor with hand inflator, three-way stopcock and colour-coded pressure indicator

Dispo Infusor

Size	With female Luer Lock	With hand inflator and three-way stopcock	With hand inflator, three-way stopcock and colour-coded pressure indicator	Box
500 ml	REF 56-01-050	REF 56-03-050	REF 56-02-050	5
1000 ml	REF 56-01-100	REF 56-03-100	REF 56-02-100	5
3000 ml	REF 56-01-300	–	REF 56-02-300	1

Infusor 500 ml / 1000 ml / 3000 ml

Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes

Features

Completely transparent

- Allows clear view on the fluid bag and the fluid level
- 360° visible

Wrap-Around Design with hook

- Fluid bag can be replaced quickly and easily
- Consistent pressure on the complete fluid bag

Manometer

- Precise pressure range from 0–300 mmHg

Hand Inflator

- Large volume allows quick inflation

Three-Way stopcock

- For pressure control
- For quick deflation

Reusable

- Easy to clean (PU material)



Infusor with female Luer Lock

Infusor with hand inflator, three-way stopcock and colour-coded pressure indicator

Infusor with hand inflator, three-way stopcock and manometer

Infusor

Size	With female Luer Lock	With hand inflator, three-way stopcock and colour-coded pressure indicator	With hand inflator, three-way stopcock and manometer	Box
500 ml	–	REF 57-02-050	REF 57-08-050	1
1000 ml	–	REF 57-02-100	REF 57-08-100	1
3000 ml	REF 57-01-300	REF 57-02-300	REF 57-08-300	1

Pressure Infusion Cuff 500 ml / 1000 ml

Pressure Infusion Cuff for rapid pressure infusion and irrigation purposes

Features

Transparent Netting

- Allows clear view on the fluid bag and the fluid level

Pocket Design with Double Hook

- Fluid bag can be replaced quickly and easily

Manometer

- Precise pressure range from 0–300 mmHg

Hand Inflator

- Large volume allows quick inflation

Three-Way stopcock

- For pressure control
- For quick deflation

Reusable



Pressure Infusion Cuff with hand inflator, three-way stopcock and colour-coded pressure indicator

Pressure Infusion Cuff with hand inflator, three-way stopcock and manometer

Pressure Infusion Cuff

Size	With hand inflator, three-way stopcock and colour-coded pressure indicator	With hand inflator, three-way stopcock and manometer	Box
500 ml	REF 58-02-050	REF 58-08-050	1
1000 ml	REF 58-02-100	REF 58-08-100	1

Pressure Infusion Cuff 3000 / 5000 ml

Pressure Infusion Cuff for irrigation purposes during minimal invasive surgery

Features

Transparent Netting

- Allows clear view on the fluid bag and the fluid level

Design

- Pocket design with zip

Manometer

- Precise pressure range from 0–300 mmHg

Hand Inflator

- Large volume allows quick inflation
- With release valve for pressure control

Reusable

- Bladder is replaceable



Pressure Infusion Cuff with female Luer Lock



Pressure Infusion Cuff with hand inflator and manometer

Pressure Infusion Cuff

Size	With female Luer Lock	With hand inflator and manometer	Box
3000 ml / 5000 ml	REF 52-01-300*	REF 52-08-300*	1

*Contains or presence of natural rubber latex.

Manometer

Pressure proved, for fixation to Pressure Infusion Cuffs



Manometer 300mmHg

REF 20-17-301

Box

1

Hand Inflator



Hand Inflator

With colour-coded pressure indicator and three-way stopcock

REF 51-02-666

Box

1



Hand Inflator

With release valve

REF 51-10-666

Box

1



Hand Inflator

With three-way stopcock

REF 51-03-666

Box

1

Automatic Pressure Infusors

For the use with Pressure Infusion Cuffs for automatic irrigation purposes during minimal invasive surgery

Features

Precision Pressure Regulator

- Very precise pressure regulation
- Maintains selected pressure and guarantees automatic compensation
- Limits pressure to max. 300 mmHg

Protection Cover

- Protects the device from liquid and damage

Coil Extension Hose

- Firmly connected to the device
- With male Luer Lock connector

Universal Clamp

- Horizontal or vertical position for attachment to drip stand or rail

Pressure Infusor 100

For compressed air, with single channel, complete with pressure hose (4 m)

Version	With pressure hose (4 m) and pressure connection RECTUS	With pressure connection NIST	Box
With universal clamp	REF 52-13-100	REF 52-13-100NIST	1
With stand and basket	REF 52-15-100	REF 52-15-100NIST	1



Pressure Infusor 200

For compressed air, with switch for the connection of two Pressure Infusion Cuffs, complete with pressure hose (4 m)

Version	With pressure hose (4 m) and pressure connection RECTUS	With pressure connection NIST	Box
With universal clamp	REF 52-13-200	REF 52-13-200NIST	1
With stand and basket	REF 52-15-200	REF 52-15-200NIST	1

